

Product Name	Cat. No.	Pack Size
--------------	----------	-----------

<b>Super Rev. Transcriptase MuLV KIT (Super RT 50 Rxns KIT)</b>	<b># BB-E0043</b>	<b>10,000 U (200 U/μl)</b>
---	-------------------	----------------------------


**Description:** Murine Leukemia Virus Reverse Transcriptase (MuLV RT) is an RNA-dependent DNA polymerase that can be used in cDNA synthesis with long messenger RNA templates (>5Kbp). The enzyme is isolated from *E.coli* expressing modified *pol* gene of MuLV on a plasmid. ***This modified version of RT is thermotolerant (working temp 42 – 52 °C using ng level of total RNA)***. The RNase H activity of MuLV RT is weaker than the commonly used AMV- Reverse Transcriptase.

**Reagents supplied with the KIT**

1. Super RT (200 U/μl), 100 μl
2. 5 x RT Buffer (250 mM Tris-HCl, pH 8.4; 375 mM KCl; 15 mM MgCl<sub>2</sub>), 400 μl
3. dNTP mix (10mM), 50 μl
4. Oligo (dT) mix (500 μg/ml), 50 μl
5. Random Hexamer (50 ng/ul), 50 ul
6. 0.1 M DTT, 100 μl
7. Murine RNase Inhibitor (1 U/μl), 50 μl
8. DEPC H<sub>2</sub>O, 600 μl

**Storage Buffer:** 20mM Tris-HCl (pH 7.5), 100mM NaCl, 0.1mM EDTA, 1mM DTT, 50% glycerol (v/v).

**Storage Instruction:** - 20°C



Reverse Transcription done with the 200 U of Super-RT (MuLV) using above mentioned amount of total RNA in a 20 μl reaction volume. 2 μl of RT product used for the PCR of 1 Kbp amplicon of GAPDH.



[www.nextgenbioproducts.com](http://www.nextgenbioproducts.com)

For more information: [contact@nextgenbioproducts.com](mailto:contact@nextgenbioproducts.com)

### **First-Strand cDNA Synthesis Using Super RT Enzyme**

- Normally, a 20  $\mu$ l Super RT reaction volume uses 5ng/10ng of Total RNA or (1ng to 2 ng of mRNA) for reaction at temperature 42°C. However, depending on the nature of reaction the input template total RNA or mRNA may be enhanced and reaction temperature may also be chosen at a workable temperature between 42°C to 52°C

- Add the following components to a nuclease-free microcentrifuge tube

<b>Use</b> Oligo (dT) <sub>12-18</sub> (500 $\mu$ g/ml) <b>or,</b> Blend of random Hexamer (50 ng / $\mu$ l ) & Oligo(dT) <sub>12-18</sub> (3:1) v/v <b>or,</b> 2 pmole gene-specific primer (GSP)	1 $\mu$ l
5 ng to 10 ng total RNA OR 1 ng to 2 ng of mRNA OR more total RNA / mRNA	X $\mu$ l
1 $\mu$ l dNTP Mix (10mM each)	1 $\mu$ l
Sterile, distilled water	Up to 12 $\mu$ l

- Heat mixture to 65°C for 5 min. & quick chill on ice. Collect the contents of the tube by brief centrifugation and add (preferably make a master mix and add accordingly):

5 X First-Strand Buffer	4 $\mu$ l
0.1 M DTT	2 $\mu$ l
RNase Inhibitor	1 $\mu$ l

- Mix contents of the tube gently. Keep it at room temp for 5 min.
- Add 1  $\mu$ l of Super RT (200U/ $\mu$ l) and mix by pipetting gently (softly) up and down
- Incubate at **42°C - 52°C** for 50 min (depending on your template)
- Inactivate the reaction by heating at 70°C for 15 min.

-----X-----