


Product Name	Cat. No.	Pack Size
DNA ladder LMW <sup>S</sup>	# BB-D60	50 µg (1ml)
	# BB-D61	2 X 50 µg (2ml)

**Concentration: 0.05 µg/µl**

**Storage: -20°C**

This ready to load low molecular weight DNA ladder contains seven bands (1412bp, 1000bp, 517bp, 396bp, 214bp, 125bp, 75/65bp) and Xylene Cyanol as tracking dye.

**Note:** 5µl (i.e. 500 ng) loading is recommended.



LMW<sup>S</sup> ladder visualized by EtBr staining on a 1 % TAE agarose gel. Mass value 500 ng / lane.